

Extended Hydration – Multipurpose Disinfection Solution (MPDS)

Reclens Extended Hydration – Multi-Purpose Disinfection Solution (MPDS) is formulated with polyhexanide, a broad-spectrum antimicrobial agent effective against a wide range of microorganisms, to ensure thorough disinfection of the lens surface. Combined with EDTA, which helps loosen deposits and build-up, it ensures a deeper clean. This solution supports daily hygiene, keeping lenses clean, fresh, and comfortable for all-day wear.



Cleaning and Disinfection
Maintains lens hygiene with antimicrobial properties that help reduce the build-up of microbes



Extended Hydration
Helps keep lenses moist for longer to address dryness and discomfort, and support comfortable lens wear throughout the day



6-in-1 Formula
New MPDS 6-in-1 formula for disinfecting, cleaning, protein removal, storing, wetting, and rinsing



Suitable for: - all types of soft contact lenses, including new generation silicone hydrogel lenses

REF	Size	Sterility	Pcs/Unit	Units/Ctn
13317013	2 x 500mL	Sterile	6	1
13317029	2 x 100mL	Sterile	10	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

