

## Lumbar Support Packs

Sage Lumbar Support Packs are sterile packs designed to provide healthcare professionals with essential tools and materials needed to support lumbar procedures efficiently and safely. Convenient and easy to use, Sage Lumbar Support Packs come in a range of three pack types offering a specialised set of components to assist with lumbar procedures for adult, infant and paediatric patients. Note: These packs do not include a spinal needle, manometer or 3-way stopcock.



### Easy Open Pack

Quick accessibility and clear labelling of contents



### Ideal for Infection Control

Packs provide an efficient on-demand sterile set up



### Various Component Options

Packs are offered in a range of configurations suitable for adult, infant and paediatric patients



11015050: Adult

### Suitable for:

- ICU
- Emergency departments
- Healthcare settings where lumbar puncture procedures are performed

Reorder details 

## Lumbar Support Packs

REF	11015050 Adult	11015051 Infant	11015052 Paediatric
<b>Sterility</b>	Sterile	Sterile	Sterile
<b>Pcs/Unit</b>	10	10	10
<b>Units/Ctn</b>	1	1	1
<b>Contains</b>			
Large Plain Tray 250mm x 200mm x 25mm	1	1	1
Syringe 3ml L/Slip CN	1	1	1
Syringe 5ml L/Lock CN	1		
Blunt Draw Up Needle ½ inch	1		
Hypo Needle 18G x 1.5" (38mm)	1		
Hypo Needle 23G x 1.25" (32mm)	1		
Hypo Needle 25G x 5/8" (16mm)	1	1	1
Gallipot 60ml	2	2	2
Non-Woven Bobs	10	10	10
Non-Woven Swabs 5cm x 5cm 4ply	10	5	5
Specitube 8ml #1	1	1	1
Specitube 8ml #2	1	1	1
Specitube 6ml #3		1	1
Specitube 8ml #3	1		
Specitube 6ml #4			1
Specitube 8ml #4	1		
Specitube 6ml #5	1		
Petri Dish	0.5	0.5	0.5
Forceps Green	2	2	2
Alcohol Swab	3		3
Strip Bandage	1	1	1
Towel Paper 42cm x 40cm 4ply	1	1	1
Plastic Drape 75cm x 100cm D/S Tape (Fen 12.5cm x 12.5cm)	1		1
Drape Blue Laminated 60cm x 30cm (Fen 4cm x 8cm)		1	
Blue Laminated Field 75cm x 75cm	1	1	1

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

