

Asguard[®]

SILICONE BORDER



Asguard Silicone Border


Asguard Silicone Border is an advanced 5 layer, all in one foam dressing designed to deliver absorption, comfort and protection across a broad range of wounds. Its gentle silicone adhesive safeguards fragile skin while preventing adherence to the wound bed, ensuring low trauma dressing changes. The highly absorbent foam core and breathable backing maintain an optimal moist wound environment, while the conformable construction adapts easily to anatomical contours and joints. Enhanced with our Omnidirectional (OmniFlex) Incision Design, the dressing provides 360° stretch for exceptional flexibility, secure wear on curved or mobile areas, and reduced risk of early edge lift; supporting reliable, extended wear without compromising comfort.



Omnidirectional Flex Technology How does it work?

The Asguard Silicone Borders main feature is the '**Omnidirectional Flex Technology**'. So how does this Tech work?

The Omnidirectional Flex Technology incorporates precisely engineered multidirectional incisions within the foam structure to enable exceptional stretch, flex, and anatomical conformability. These incisions also optimise exudate transmission pathways, providing a clearer and more reliable strikethrough signal to support timely dressing changes.



What shapes does the Silicone Border come in?



Rectangle Dressing



Square Dressing



Heel Dressing



Sacral Dressing

Extended, Secure, Gentle

Asguard Silicone Border


Engineered with an extended silicone border, the dressing delivers superior stay on performance by maximising skin contact and boosting overall adhesion. The broader silicone footprint significantly reduces edge roll and lift, helping the dressing remain secure for longer, even on high movement areas.

Compared to traditional adhesives or contact layers, Asguard Silicone high-adhesive technology offers a balance of secure fixation, gentle removal, and clinical versatility. It is especially beneficial for:

- Chronic wounds
- Fragile or sensitive skin
- Pediatric and geriatric patients
- Post-surgical sites



Our layer breakdown for Silicone Border



Omnidirectional Flex Technology

The incisions within the foam structure allows for exceptional stretch, flex, and conformability



Asguard Silicone Adhesive

High adhesion featuring a larger Silicone Border for improved securement and extended wear time



Absorbent Foam Core

Manages exudate efficiently



Water & Shower Resistant Film

Designed to reduce dressing changes and maintain integrity during routine hygiene and mobility

* The specified goods should only be used under the supervision of a healthcare professional.

* If you are concerned about your wound, consult a healthcare professional



REF	Dressing Size	Padding Size	Equivalent To	Type	Sterility	Pcs/Box	Boxes/Ctn
AWS575	8.9cm x 8.9cm	4.5cm x 4.5cm	7.5cm x 7.5cm	Square	Sterile	10	20
AWS50512	8.2cm x 15.3cm	3cm x 8.7cm	7.5cm x 12.5cm	Rectangle	Sterile	10	20
AWS510	11.55cm x 11.55cm	6.3cm x 6.3cm	10cm x 10cm	Square	Sterile	10	20
AWS512	14.2cm x 14.2cm	9cm x 9cm	12.5cm x 12.5cm	Square	Sterile	10	20
AWS515	16.8cm x 16.8cm	11.6cm x 11.6cm	15cm x 15cm	Square	Sterile	5	40
AWS518S	18cm x 18cm	19cm x 20cm	-	Sacral	Sterile	5	40
AWS523S	23cm x 23cm	12.5cm x 14cm	-	Sacral	Sterile	5	40
AWS523H	23.9cm x 23.9cm	17cm x 18cm	-	Heel	Sterile	5	40