

## Silicone Border

Asguard Silicone Border is an all-in-one 5-layer dressing designed to deliver absorption, comfort and protection across a broad range of wounds. Its soft silicone adhesive layer gently adheres to the skin without sticking to the wound bed, minimising pain and trauma during dressing changes. The dressing features a highly absorbent foam pad that manages exudate effectively, while the breathable backing film maintains an optimal moist wound environment.



### Extended Silicone Border

Expanded silicone border around the absorbent layers reduces the risk of rolling or detachment during use



### Omnidirectional Flex Technology\*

The incisions within the foam structure allows for exceptional stretch, flex, and conformability  
\*Incisions not available on the sacral or heel dressings



### Water & Shower Resistant Film

Designed to reduce dressing changes and maintain integrity during routine hygiene and mobility activities



### 5-Layer Design

Advanced 5-layer dressing design for effective wound management



- Suitable for:**
- Superficial wounds with light to moderate exudate
  - Arterial, diabetic and venous leg ulcers
  - Skin tears, pressure injuries, partial thickness burns and surgical wounds

REF	Type	Dressing Size	Pad Size	Equivalent To	Sterility	Pcs/Box	Boxes/Ctn
<b>AWS575</b>	Square	8.9cm x 8.9cm	4.5cm x 4.5cm	7.5cm x 7.5cm	Sterile	10	20
<b>AWS510</b>	Square	11.55cm x 11.55cm	6.3cm x 6.3cm	10cm x 10cm	Sterile	10	20
<b>AWS512</b>	Square	14.2cm x 14.2cm	9cm x 9cm	12.5cm x 12.5cm	Sterile	10	20
<b>AWS515</b>	Square	16.8cm x 16.8cm	11.6cm x 11.6cm	15cm x 15cm	Sterile	5	40
<b>AWS50512</b>	Rectangle	8.2cm x 15.3cm	3cm x 8.7cm	7.5cm x 12.5cm	Sterile	10	20
<b>AWS518S</b>	Sacral	18cm x 18cm	19cm x 20cm	-	Sterile	5	40
<b>AWS523S</b>	Sacral	23cm x 23cm	12.5cm x 14cm	-	Sterile	5	40
<b>AWS523H</b>	Heel	23.9cm x 23.9cm	17cm x 18cm	-	Sterile	5	40

Asguard® Silicone Border should only be used under the supervision of a healthcare professional. If you are concerned about your wound, consult a healthcare professional.

POUCH - PAPER



POUCH - FILM



BOX



**Sentry Medical Pty Ltd**

1300 995 999 | sales@sentrymedical.com.au  
www.sentrymedical.com.au

**Sentry**  
MEDICAL