

Compostable Apron

Owear Compostable Aprons are single use and compostable. The Owear Compostable Apron provides fluid barrier protection of a typical disposable apron with less environmental impact by conforming to DIN EN 13432 and ASTM D 6400 compostability standards.



Industrially Compostable

*Meets DIN EN 13432 and ASTM D 6400 compostable product test standards



Waterproof

Protection against fluids, spills, and splashes



Easy Tie

Neck loop and back ties for ease of donning and doffing

Suitable for:

- Protection against fluid splashes
- Protection against dirt and dust
- Protection against general grease and grime



REF	Size	Colour	Pcs/Unit	Units/Ctn
AP001C	61cm x 122cm	White	100	10
AP002C	81cm x 132cm	White	100	6
AP003C	81cm x 142cm	White	100	6
AP004C	81cm x 160cm	White	100	6

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

