

ASTM Level 2 Ultra Soft Earloop Masks

PRIMED® have engineered a range of ASTM Level 2 Masks with unique Ultra Soft Earloops for added comfort, making them softer and thicker. These ASTM Level 2 masks provide a moderate level of protection during medical procedures, making them a go-to mask for basic protection during wound care, line changes, dental work and procedures with moderate risk of fluid splash.



BFE (%): ≥ 98
 PFE (%): ≥ 85
 ΔP (mmH₂O/cm²): <6.0
 Fluid Resistance: 120mm Hg
 Flame Spread: Class 1

Variation in lot results exist.
 All results must meet or exceed
 ASTM F2100 Requirements.



Ultra Soft Earloop

Experience ultimate comfort with our premium ear loops



Ultrasonically Welded Seams

Ensures mask layers remain secure and rest comfortably against the clinician's skin



Wide Nosepiece

Constructed with high-quality aluminium to keep the mask in place across the nose and cheekbones



Soft Mask

Enhances comfort with its soft, hypoallergenic inner layer, ideal for those with sensitivities to other masks*

*Applies to PG4-1312, PG4-1342, PG4-1352.



REF	Feature	Hypoallergenic	Colour	Pcs/Unit	Units/Ctn
PG4-1832	Ultra Soft Earloops	-	Blue	50	10
PG4-1872	Anti Fog Foam Mask, Ultra Soft Earloops	-	Blue	50	10
PG4-1342	Anti Fog Foam Soft Mask, Ultra Soft Earloops	Yes	White	50	10
PM4-314	Anti Fog Strip Mask, Ultra Soft Earloops	-	Blue	50	10
PG4-1312	Soft Mask, Ultra Soft Earloops	Yes	White	50	10
PG4-1352	Soft Mask, Ultra Soft Earloops	Yes	Blue	50	10

At Sentry, we are dedicated to upholding sustainable principles to safeguard our planet, people, and future, and we strive to consistently act responsibly in doing so.

Sentry Medical Pty Ltd

1300 995 999 | sales@sentrymedical.com.au
 www.sentrymedical.com.au

