Asguard®

Hydroporous Foam

Asguard® Hydroporous Foam - Non-Adhesive is an ideal dressing to manage wounds with different needs and moisture levels. The fine pores retain moisture at the wound surface, while excess exudate travels into upper pores of the foam to be eliminated as moisture vapour, increasing dressing wear time.

Moisture Vapour Permeable

Allows moisture vapour transmission while maintaining a moist wound environment, helping to manage exudate and reduce the risk of peri-wound maceration

Waterproof Backing

Waterproof outer film provides a barrier to water and bacteria while remaining breathable

Contours to Body

Soft, conformable foam adapts to body contours, helping maintain close contact with the wound surface for improved comfort and effective exudate management



Suitable for:

- Primary and/or secondary dressing for the management of wounds
- Arterial, diabetic and venous leg ulcers
- Pressure injuries, partial thickness burns and surgical incisions

REF	Description	Size	Pcs/Unit	Units/Ctn
AWF250	Hydroporous Foam - Non-Adhesive	5cm x 5cm	10	10
AWF251	Hydroporous Foam - Non-Adhesive	6cm x 8.5cm	5	10
AWF252	Hydroporous Foam - Non-Adhesive	10cm x 10cm	10	10
VARI-SITE				
AWF260	Hydroporous Foam - Non-Adhesive	12cm x 19cm	5	10

Asguard® Hydroporous Foam - Non-Adhesive should only be used under the supervision of a healthcare professional. If you are concerned about your wound, consult a healthcare professional.

POUCH - PAPER	POUCH - FILM	BOX
RECYCLE ME	GENERAL WASTE	RECYCLE ME

Sentry Medical Pty Ltd
1300 995 999 | sales@sentrymedical.com.au
www.sentrymedical.com.au

Sentry MEDICAL

SMHO-MKT-FLY-0100-REV4_10_26 Seek advice from a health professional before use. All product information is correct at date of publication.